



# HICHROM

Chromatography Columns and Supplies

## APPENDICES HPLC Calculations


Catalogue 9

### Hichrom Limited

1 The Markham Centre, Station Road  
Theale, Reading, Berks, RG7 4PE, UK

Tel: +44 (0)118 930 3660 Fax: +44 (0)118 932 3484

Email: [sales@hichrom.co.uk](mailto:sales@hichrom.co.uk) [www.hichrom.co.uk](http://www.hichrom.co.uk)



where:

- $t_0$  = Retention time of unretained peak
- $t_{R1}$  = Retention time of component 1
- $t_{R2}$  = Retention time of component 2
- $w_{0.5}$  = Peak width at 50% peak height
- $w_{0.1}$  = Peak width at 10% peak height
- $w_{0.05}$  = Peak width at 5% peak height
- $w$  = Peak width at base
- $A$  = Peak front distance at 10% of peak height
- $B$  = Peak tail distance at 10% of peak height
- $C$  = Peak front distance at 5% of peak height
- $D$  = Peak tail distance at 5% of peak height

The following equations are based on the terms defined above.

Adjusted (net) Retention Time ( $t'_R$ )  $t'_{R1} = t_{R1} - t_0$

### Selectivity

Retention Factor ( $k$ )  $k = \frac{t_R - t_0}{t_0} = \frac{t'_R}{t_0}$

Separation Factor ( $\alpha$ )  $\alpha = \frac{k_2}{k_1}$

### Performance

Column Efficiency ( $N$ )  
(plate number or number of theoretical plates)  $N_{0.5} = 5.54 \left( \frac{t_R}{w_{0.5}} \right)^2$  or  $N_{0.1} = 18.55 \left( \frac{t_R}{w_{0.1}} \right)^2$

Height Equivalent to a Theoretical Plate (HETP)  $H = \frac{L}{N}$  where  $L$  = length of column

Reduced Plate Height ( $h$ )  $h = \frac{H}{d_p}$  where  $d_p$  = average particle diameter

### Peak Shape

Asymmetry Factors (see p.14)  $As_1 = \frac{N_{0.1}}{N_{0.5}}$       $As_2 = \frac{B}{A}$

USP Tailing Factor  $T = \frac{W_{0.05}}{2C}$

Resolution  $R_s = 2 \left( \frac{t_{R2} - t_{R1}}{w_{1+2}} \right)$

or

$R_s = \frac{\sqrt{N}}{4} \left( \frac{\alpha - 1}{\alpha} \right) \left( \frac{k}{k+1} \right)$  where  $k$  is the average retention factor for peaks 1 and 2