



HICHROM

Chromatography Columns and Supplies

**LC COLUMNS
SIELC
SHARC**

Catalogue 9

Hichrom Limited

**1 The Markham Centre, Station Road
Theale, Reading, Berks, RG7 4PE, UK**

Tel: +44 (0)118 930 3660 Fax: +44 (0)118 932 3484

Email: sales@hichrom.co.uk www.hichrom.co.uk

SHARC™

SHARC™ (Specific Hydrogen-bond Adsorption Resolution Chromatography) columns are the latest columns developed by SIELC Technologies. They are the first commercially available columns which provide separations based primarily on hydrogen bonding. SHARC 1 achieves separation based on the analyte's ability to act as a hydrogen atom donor or acceptor.

Operating Conditions

Solvents used for SHARC separations are acetonitrile (MeCN) as the weak solvent and methanol (MeOH) as the strong solvent. Acetonitrile alone has a low level of hydrogen bonding with the SHARC stationary phase, whereas methanol interacts strongly, reducing retention of

analytes capable of hydrogen interactions (Figure 16). By altering the ratio of acetonitrile/methanol, the optimum retention profile can be obtained for many types of molecules.

Analytes can retain on the stationary phase by more than one hydrogen bond and act as a donor or acceptor. The SHARC 1 column is a hydrogen atom acceptor type stationary phase, showing increased retention towards molecules with higher numbers of polar X-H bonds such as alcohols, amines, acids, amides, phenols etc.




Figure 16. SHARC technology

Benefits of SHARC Columns

- 1) Acetonitrile-methanol mixtures have lower viscosity than aqueous solvent mixtures, enabling the use of smaller particles and faster flow rates, for fast analyses.
- 2) Methanol is one of the most universal solvents for organic compounds.
- 3) Acetonitrile-methanol mixtures have a low boiling point and are much easier to evaporate than water, a benefit for preparative separations.
- 4) Acetonitrile-methanol mixtures are MS friendly, enabling mass directed preparative strategies.
- 5) A wide range of compounds with functional groups containing oxygen and nitrogen can be retained and separated.

Figure 17 shows the separation of a mixture of four neurotransmitters, based on their ability to form hydrogen bonds with the SHARC 1 stationary phase. Elution order corresponds to the number and strength of interaction points.




Figure 17. Effect of eluent composition on separation on SHARC 1

Ordering Information – SHARC 1 (5µm, 100Å)

Column i.d. ¹ (mm)	Column Length ¹ (mm)				Guards ² (2/pk)
	50	100	150	250	
2.1	SH1-21.050.0510	SH1-21.100.0510	SH1-21.150.0510	SH1-21.250.0510	SH1-21.G.0510
3.2	SH1-32.050.0510	SH1-32.100.0510	SH1-32.150.0510	SH1-32.250.0510	SH1-32.G.0510
4.6	SH1-46.050.0510	SH1-46.100.0510	SH1-46.150.0510	SH1-46.250.0510	SH1-46.G.0510
10	SH1-100.050.0510	SH1-100.100.0510	SH1-100.150.0510	SH1-100.250.0510	-
22	SH1-220.050.0510	SH1-220.100.0510	SH1-220.150.0510	SH1-220.250.0510	-
22 (10µm)	SH1-220.050.1010	SH1-220.100.1010	SH1-220.150.1010	SH1-220.250.1010	-

¹ Other dimensions available, including capillary

² Direct connect - no holder required