

HPLC Column for hydrophilic interaction chromatography

COSMOSIL HILIC Packed Column


COSMOSIL HILIC is a new column for hydrophilic interaction chromatography with Triazole bonded silica packing material. The hydrophilic interaction chromatography is a variation of normal phase chromatography where a polar stationary phase is used with a mobile phase which contains a high concentration of organic solvent and a low concentration of aqueous eluent. The main retention mechanism is the partitioning of the polar analytes between the polar stationary and the non-polar mobile phase. As it is also called "aqueous normal phase", the elution order is similar to that of normal phase and the sample elution is in the order of increasing hydrophilicity. Without using ion-pair reagent COSMOSIL HILIC retains highly polar analytes that would not be retained in reversed phase chromatography. It also shows a weak anion-exchange mechanism with the positively charged stationary phase, thus acidic compound is strongly retained.

Features of COSMOSIL HILIC

- Retains highly polar compounds that would not be retained in C18 columns.
- Ion-pair reagent is not required
- Anion-exchange results in strong retention for acidic compounds.
- Enhanced sensitivity in LC/MS


• Comparison with C18 Columns

COSMOSIL HILIC can separate Glycine and Glycylglycine without ion-pair reagent. Although C18 column can separate them with ion-pair reagents, there are some disadvantages such as column equilibration, preparation of mobile phase and column deterioration.





Material characteristics

Packing material	COSMOSIL Cholester
Silica gel	high purity spherical porous silica
Average particle size	5µm
Average pore size	12 nm
Specific surface area	approx. 300 ² /g
Stationary phase	Triazole
Main interaction	Hydrophilic interaction



Different interactions


The retention Mechanism of COSMOSIL HILIC is the combination of hydrophilic interaction and anion exchange, and the retention can be controlled by changing the mobile phase.



Analysis data


BCAA (amino acid branched-chain)

Column size 4.6mm I.D. x 250mm
 Mobile phase Acetonitrile:10mmol/l CH₃COONH₄=85:15
 Flow rate 1.0ml/min
 Temperature 30°C
 Detection ELSD (Atten=8, Gain=6)
 Sample 1. Leucine (1.0mg/ml)
 2. Isoleucine (1.0mg/ml)
 3. Valine (1.0mg/ml)
 Injection volume 3.0µl




Juice components

Column size 4.6mm I.D. x 250mm
 Mobile phase Acetonitrile:10mmol/l Phosphate Buffer (pH7.0)=50:50
 Flow rate 1.0ml/min
 Temperature 30°C
 Detection UV 210nm
 Sample 1. Ascorbic Acid (0.5mg/ml)
 2. Malic Acid (1.0mg/ml)
 Injection volume 3.0µl




Phosphorylated saccharide

Column size 4.6mm I.D. x 250mm
 Mobile phase Acetonitrile:20mmol/l Phosphate Buffer (pH7.0)=60:40
 Flow rate 30°C
 Temperature RI
 Detection 1. Glucose (2mg/ml)
 Sample 2. α-D-Glucose-1-Phosphate (10mg/ml)
 3. D-Glucose-6-Phosphate (10mg/ml)
 Injection volume 5.0µl




Water-soluble vitamin

Column size 4.6mm I.D. x 250mm
 Mobile phase Acetonitrile:100mmol/l CH₃COONH₄=80:20
 Flow rate 1.0ml/min
 Temperature 30°C
 Detection UV 220nm
 Sample 1. Nicotinamide (0.25mg/ml)
 2. Pyridoxine (Vitamin B₆) (0.50mg/ml)
 3. Vitamin B₂ (Riboflavin) (0.50mg/ml)
 4. Nicotinic Acid (0.25mg/ml)
 5. D-Pantothenic Acid (6.25mg/ml)
 6. L-(+)-Ascorbic Acid (1.75mg/ml)
 Injection volume 0.5µl




Hydrophilic drug

Column size 4.6mm I.D. x 250mm
 Mobile phase Acetonitrile:10mmol/l CH₃COONH₄=90:10
 Flow rate 1.0ml/min
 Temperature 30°C
 Detection UV 254nm
 Sample Famotidine (0.25mg/ml)
 Injection volume 2.0µl




Natural product

Column size 4.6mm I.D. x 250mm
 Mobile phase Acetonitrile:10mmol/l CH₃COONH₄=80:20
 Flow rate 1.0ml/min
 Temperature 30°C
 Detection UV 220nm
 Sample L-(+)-Ergothioneine (1.0mg/ml)
 Injection volume 3.0µl




GABA(4-Amino-n-butyric Acid)

Column size 4.6mm I.D. x 250mm
 Mobile phase Acetonitrile:10mmol/l CH₃COONH₄=75:25
 Flow rate 1.0ml/min
 Temperature 30°C
 Detection UV 210nm
 Sample 1. GABA(4-Amino-n-butyric Acid) (10mg/ml)
 Injection volume 3.0µl




Cyanoacetic Acid

Column size 4.6mm I.D. x 250mm
 Mobile phase Acetonitrile:10mmol/l CH₃COONH₄=50:50
 Flow rate 1.0ml/min
 Temperature 30°C
 Detection ELSD
 Sample 1. Cyanoacetic Acid (10mg/ml)
 Injection volume 0.5µl



Combination with C18 columns




Selection of mobile phase

Following are the recommended mobile phases for different samples.

- Neutral compounds Acetonitrile / Water = 90 / 10
- Basic compounds Acetonitrile / 10mmol/l CH₃COONH₄ = 90 / 10
- Amphoteric compounds Acetonitrile / 10mmol/l CH₃COONH₄ = 70 / 30
- Acidic compounds Acetonitrile / 10mmol/l CH₃COONH₄ = 50 / 50
 (not eluted) → Acetonitrile / 10mmol/l Phosphatase buffer (pH7.0) = 50 / 50

COSMOSIL Chromatogram Index

COSMOSIL HILIC Chromatogram Index, which includes 154 chromatograms using COSMOSIL HILIC, is now available. This index is useful for optimizing analytical conditions for hydrophilic interaction chromatography.



Ordering information

Product name	Size	Code No.	Product name	Size	Code No.
COSMOSIL HILIC	2.0mm I.D. x 50mm	07052-91	COSMOSIL HILIC	4.6mm I.D. x 10mm	07055-61
Packed Column	2.0mm I.D. x 150mm	07054-71	Guard Column	10.0mm I.D. x 20mm	07058-31
	4.6mm I.D. x 150mm	07056-51			
	4.6mm I.D. x 250mm	07057-41			
	10.0mm I.D. x 250mm	07059-21			
	20.0mm I.D. x 250mm	07060-81			

Other size may be available. Please enquire!!

For research use only, not intended for diagnostic or drug use.



Hichrom Limited

1 The Markham Centre, Station Road
 Theale, Reading, Berks, RG7 4PE
 Tel: 0118 930 3660
 Fax: 0118 932 3484
 Web site: www.hichrom.co.uk
 Email: sales@hichrom.co.uk

NACALAI TESQUE, INC.

Nijo Karasuma, Nakagyo-ku, Kyoto 604-0855, JAPAN
 TEL: +81-75-251-1730
 FAX: +81-75-251-1763
 Web site: http://www.nacalai.com
 E-mail: info.intl@nacalai.co.jp